

# SHARPSGUARD<sup>®</sup>

## Container in use



1

Ensure container is situated at the point of use



2

Slide door to temporary closure in-between use



3

Lock container by sliding door to final closure position and complete label for disposal



If you have questions or need further information, please contact your local Area Manager or our Customer Service Team on +44 1706 754980 or [info.healthcareuk@mauserpackaging.com](mailto:info.healthcareuk@mauserpackaging.com) and visit our website [www.daniels.co.uk](http://www.daniels.co.uk).



[www.mauserpackaging.com](http://www.mauserpackaging.com)



[info@mauserpackaging.com](mailto:info@mauserpackaging.com)

REDEFINING SUSTAINABILITY

Unless otherwise stated, the information contained herein is at the specific request of the user and has been made available by Daniels Healthcare.

V01.2023