

# SHARPSGUARD®

## Container in Use



1

Ensure container is situated at the point of use



2

Slide door to temporary closure in-between use



3

Lock container by sliding door to final closure position and complete label for disposal



WILTON  
HEALTH  
CARE

**MAUSER**  
Packaging Solutions

Email: [sales.dept@wilton-hc.co.uk](mailto:sales.dept@wilton-hc.co.uk) Tel: 028 9022 4433

Safety, Service and Support - *Every Day, Affordably & Responsibly*

[www.daniels.co.uk](http://www.daniels.co.uk)

**d.**  
DANIELS®

A Mauser Packaging Solutions brand

NI Container in Use V2-2024