

# SHARPSGUARD® ECO

## Theatre+ Assembly

1

Place lid on the container with the opening at the front



Lid Opening

Label

Feet

2

Press firmly at the front and back of the container



3

Press firmly at each side of the container



4

Check the corners are in place by pressing firmly on each corner



5

Repeat step 4 on the opposite corners



If you have questions or need further information, please contact your local Area Manager or our Customer Service Team on +44 1706 754980 or [info.healthcareuk@mauserpackaging.com](mailto:info.healthcareuk@mauserpackaging.com) and visit our website [www.daniels.co.uk](http://www.daniels.co.uk).