

Category A Waste

UN3549

We have two 60L WIVA™ Containers with solid lids for use with two UN bulk approved yellow clinical waste bags which meet the requirements for UN3549 Category A Waste.

There are the:

Yellow base – Virgin material **Grey base** – Recycled material

These standard 60L WIVA™ with a solid lid are available via your normal supply routes.

INSTRUCTIONS FOR USE

- **Inner packaging** – UN bulk approved yellow clinical waste bag; bag to be closed per bag approval (usually swan necked and bag tied) and placed inside the intermediate packaging. If free liquid present, absorbent material, like gelling agent, to be placed in the bottom of inner bag
- **Intermediate packaging** - UN bulk approved yellow clinical waste bag; bag to be closed per bag approval (usually swan necked and bag tied) and placed in the outer packaging
- **Outer packaging** – 60L WIVA™ container with a solid yellow lid. If free liquid present, absorbent material, like gelling agent, to be placed in the bottom of outer container
- **Attach UN3549 label to container** overlaying UN3291 label – (Blue panel for yellow base, Green panel for grey base) – Label available to download here: daniels.co.uk/resources/brochures/

Weight limit of 15kg



If you have questions or need further information, please contact your local Area Manager or our Customer Service Team on +44 (0)1865 371841 or Orders.UK@mauserpackaging.com and visit our website www.daniels.co.uk.